EN - IMPORTANT MEDICAL INFORMATION - RI.KNEE ADJUSTABLE TIBIA CUT GUIDE

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MATERIALS AND INTENDED USE

RI.KNEE Adjustable Tibia Cut Guide is supplied by Enztec Limited (Enztec) and is a surgical instrument intended for transient use during orthopaedic surgery. It is intended to be used in conjunction with the CORI surgical system during total knee arthroplasty. The device is manufactured from medical grade metals and plastics and is supplied non-sterile. Prior to use, the operating surgeon shall have given careful consideration to all aspects of the surgical intervention as well as to the limitations of the device.

For REAL INTELLIGENCE CORI instructions for use, refer to the REAL INTELLIGENCE CORI for Knee Arthroplasty User Manual (CORI User Manual). The instruments are compatible with all Smith+Nephew tibia baseplates supported by CORI, except bi-cruciate retaining baseplates.

INDICATIONS FOR USE

RI.KNEE Adjustable Tibia Cut Guide is intended to be used in conjunction with the CORI surgical system during total knee arthroplasty.

EXAMINATION PRIOR TO USE

The device is supplied non-sterile, in device containers or individually packaged. Device packaging must be intact when received and removed prior to sterilisation.

The device must be carefully and completely examined for wear or damage by doctors and staff in operating centres prior to surgery. The examination shall include a visual and functional inspection of the working surfaces and must verify the cleanliness of the device and the absence of wear and damage. Evidence of damage or wear on a device may include but is not limited to corrosion, discoloration, excessive scratches, flaking, distortion, wear, cracks, and loosening of components. Improperly functioning devices, devices with unrecognisable markings, and damaged or excessively worn devices should not be used.

Like any precision surgical instrument, the device must undergo regular checks by authorised personnel to ensure that the device remains in good condition and continues to act as intended. End of useful instrument life is generally determined by wear or damage in surgical use.

ASSEMBLY**Components**

DESCRIPTION	CATALOGUE ITEM
RI.KNEE Tibia Adjustment	4240-7500
RI.KNEE Tibia Fixation Base Left	4240-7501
RI.KNEE Tibia Fixation Base Right	4240-7502
RI.KNEE Tibia Cut Guide Left	4240-7503
RI.KNEE Tibia Cut Guide Right	4240-7504
RI.KNEE Tibia Tracker Adapter	4240-7505

Select the correct Tibia Fixation Base and Tibia Cut Guide depending on the operative side. Each side-specific piece of hardware is designated LEFT or RIGHT. Rotate the right dial (FLEX) counter-clockwise to the unlock symbol position (1). Insert the Tibia Cut Guide (2) then rotate the right dial (FLEX) to the 0 position (3). Insert the rack on the Tibia Adjustment into the rectangular hole in the Tibia Fixation Base (4) and rotate the screw clockwise to the 0 position (5). Prior to cleaning and sterilisation, reverse this process to fully disassemble.


Figure 1: Assembly of Tibia Adjustable Cut Guide

Unscrew the lock knob (1) on the Tracker Adapter and insert the end of the Tibia Tracker (2). Rotate the tracker to the desired position then tighten the lock knob (3). Optionally use the Checkpoint Tool to tighten. Attach Tracker Adaptor to the Fixation Base (4).


Figure 2: Attach Tracker with Quick Connect Mechanism

INSTRUCTIONS FOR USE

Follow the instructions in the CORI User Manual for beginning a CORI case, setting up the system and instruments, and preparing the patient, including making the incision. Install the Femur Tracker and attach Flat Markers to the Tibia Tracker as instructed in the CORI User Manual.

Installing Tibia Fixation Base

The Tibia Adjustment allows for $\pm 6\text{mm}$ resection, $\pm 8^\circ$ varus/valgus, and $\pm 8^\circ$ AP slope, from the zero-position of the cut slot, as described in Figure 3.


Figure 3: Adjustment Ranges

To install the Adjustable Tibia Cut Guide, ensure all dials are set to 0 to allow for the full range of adjustment. The cut slot should be aligned with the anticipated resection. Install the Tibia Fixation Base with the most lateral pin on or medial to the tibia tubercle as shown in Figure 4. Optionally, the following manual alignment instruments are compatible for use with the Adjustable Tibia Cut Guide. Refer to the associated surgical technique guide for instructions for use.

	DESCRIPTION	CATALOGUE ITEM
SYNC INSTRUMENTATION	Tibial Stylus Adjustable Footed	7401-6229
	Extramedullary Tibial Alignment Tube Macro	7401-6232
	Extramedullary Posterior Slope Guide	7401-6231
	Ankle Clamp	7401-6241


Figure 4: Cut Guide Fixation and Pin Placement

Once alignment is established, ensure the pin locations are located within the incision. Position the most lateral pin guide of the Fixation Base on the tibial tubercle medial of the subluxed patellar tendon and ensure that the most medial pin guide of the Fixation Base is within the wound. First place the most lateral pin (1), followed by the most medial pin (2) as shown in the figure above, ensuring adequate purchase in the bone. For small patients, an alternate medial pin (3) placement is available as shown in Figure 4. If further stability is desired, a third pin can be placed (4).

	DESCRIPTION	CATALOGUE ITEM
SPEED PINS FOR OPTIMAL CUTTING BLOCK FIXATION	65mm Non-Rimmed SPEED PIN	7401-3480
	80mm Non-Rimmed SPEED PIN	7401-3490

Insert the Tibia Tracker Adapter into the Tibia Fixation Base. With the Flat Markers oriented towards the Tracking Camera, secure the Tibia Tracker in the Tibia Tracker Adapter. After completing the Camera Orientation Adjustment stage, tighten the Tracker into the Adapter.

Inserting Checkpoint Pins

Insert the femur Checkpoint Pin as instructed in the CORI User Manual. Install the tibia Checkpoint Pin as follows:

- Place close to the most medial speed pin.
- Avoid interference with the final location of the cutting block.
- Consider the line of sight of the Point Probe to the Tracking Camera during checkpoint verification.


Figure 5: Placement of Checkpoint Pin

Planning and Placement of Tibial Implant

Complete the CORI registration and planning, including femur and tibia free collection, gap assessment, and femur and tibia component placement as instructed in the CORI User Manual. On the Tibia Cut Selection screen,

select Bur All or Twin Peg. Verify checkpoints as specified in the CORI User Manual. If needed, the Tibia Cut Guide can be temporarily removed to facilitate access to the nearby bone.

Upon entry into the Tibia Bone Removal screen, select Visualize Cut as described in the CORI User Manual. Attach the Plane Visualization Tool to the Point Probe and insert into the cut slot, ensuring contact with the bottom surface of the slot. Align the Plane Visualization Tool with the planned cut, using the hex screwdriver (Catalogue Item 115035) to adjust the mechanism dials that will alter the placement of the cut slot as follows:

RESECTION PLANE ADJUSTMENT	DIAL	ADJUSTMENT
	Left dial (VAR)	Varus/valgus
	Right dial (FLEX)	AP Slope
	Lower dial	Resection depth

When satisfied with the adjustment, secure the Tibia Cut Guide using speed pins. After pin installation, use the Plane Visualization Tool to confirm proper placement of the cut slot in relation to the plan.

Using a 1.35mm saw blade, perform the tibia cut through the cut slot. Use the Plane Visualization Tool to confirm the cut. If desired, refine the bone model using the Robotic Drill to determine if additional resection is needed through use of the bur.

Complete PostOp Gap Assessment using the trial components. Remove the Adjustable Tibia Cut Guide prior to preparing the keel, as the tibia pins may interfere with the tibial keel.

Once removed, fully disassemble Adjustable Tibia Cut guide prior to cleaning and sterilisation.

WARNINGS AND PRECAUTIONS

- Enztec devices must only be used by surgeons who are fully familiar with the surgical technique required and who have been trained to this end.
- The operating surgeon must take care not to exert inappropriate stress on the device and must fully comply with the operating procedure described in the Surgical Protocol.
- A full, conventional instrumentation tray for the chosen implant should be available during every use of CORI in case the surgeon must revert to a manual procedure.

Incorrect maintenance, cleaning or handling may render the device unsuitable for its intended use, cause corrosion, dismantling, distortion and/or breakage, or cause injury to the patient or operating staff. Enztec shall not be responsible in the event of a device being used which is damaged, incomplete, showing signs of excessive wear and tear, or that has been repaired or modified (either permanently or temporarily) outside the control of Enztec or its representatives.

As a result of mechanical features required, the device is made of non-implantable materials. In the event of the device breaking, no fragment must remain in the patient as this could cause post-operative complications such as allergies, infections, or complications of a biological nature associated with the release of non-implantable components, possibly requiring further intervention.

These instruments are not intended for use in a MR (Magnetic resonance) environment.

STORAGE

The device should be stored in individual packages or in containers. After use they must be stored in a clean, dry and temperate place.

DISPOSAL

Devices must be disposed of in accordance with the health care facility's procedures, ensuring protection from physical hazards such as exposed edges. Care must be taken to ensure that used devices are decontaminated following the cleaning and sterilisation instructions described within this document, or otherwise disposed of as infectious waste. Devices should be destroyed in a manner that prevents potential reuse.

COMPLAINTS

Any health professional having a complaint or grounds for dissatisfaction relating to the quality of the product, its identity, durability, reliability, safety, effectiveness, and or its performance, should notify Enztec or their representatives. Moreover, if the device has malfunctioned, or is suspected of having malfunctioned, Enztec or their representative must be advised immediately.

If an Enztec product has ever worked improperly and could have caused or contributed to a serious incident, serious injury, or death, Enztec or their representative, as well as the competent authority of the Member State in which the user and/or patient is established, must be informed as soon as possible by telephone or in writing.

For all complaints, please include the device name and catalogue number, a full description of any markings, contact name and address, and an exhaustive description of the event to help Enztec understand the causes of the complaint. Please retain the device for investigation purposes.

CLEANING AND STERILISATION


Figure 6: Tray Configuration

For safety reasons, non-sterile devices must be pre-cleaned, cleaned (manual cleaning OR automated cleaning), and sterilised before use. Moreover, for good maintenance, reusable devices must be pre-cleaned, cleaned (manual cleaning OR automated cleaning), and sterilised after surgery following the sequence of steps described hereafter:

Point of Use	- Remove gross contamination			
Transport to Processing Area	- Avoid damage - Minimise time before cleaning			
Preparation for Cleaning	- Disassemble where possible following disassembly instructions, as applicable - Instruments must be cleaned separately from instrument trays and cases			
Pre-cleaning	- Submerge in enzymatic detergent prepared according to manufacturer's recommendations - Soak for 10 minutes at 40°C (104°F) - Scrub while submerged with soft sponge and agitate - Use pipe cleaner or non-metallic brush for lumens and crevices. Actuate moving parts to loosen trapped soil - Rinse in 38°C-49°C (101°F-120°F) tap water for 1 minute - Thoroughly flush all lumens and difficult to reach areas - Actuate while rinsing			
Manual Cleaning	- Soak in ultrasonic bath in neutral pH detergent (Neutrad or acceptable alternative), prepared according to manufacturer's recommendations - Clean for 10 minutes at 40°C (104°F) - Rinse with clean, tap water, actuating moving parts while rinsing, for 1 minute - Dry thoroughly with clean, lint-free cloth			
Automated Cleaning	Phase	Recirculation Time	Water Temp	Detergent type
	Pre-wash	02:00	Cold Tap Water	N/A
	Enzyme wash	02:00	Hot Tap Water	Enzymatic Cleaner
	Wash	02:00	65.5°C (150°F)	Neutral pH Detergent
	Rinse	01:00	Hot Tap Water	N/A
	*Thermal Decontamination	01:00	90°C (194°F)	N/A
	Dry	07:00	115°C (239°F)	N/A
	*Thermal decontamination may be performed as an additional step to render the devices safe for handling			
Inspection	- Inspect for contaminants and damage - Repeat cleaning if contaminants remain - Contact Enztec or its representatives if device is damaged			

Preparation for Sterilisation	<ul style="list-style-type: none"> - Reassemble where possible following assembly instructions - Instruments with articulating surfaces must be tested for movement. A moist heat compatible, medical grade lubricant should be applied to all articulating joints prior to sterilisation (Rudolf Oil Spray RU8880-00 or acceptable alternative) - Place cleaned, dry devices into the specified locations within the cases provided, if applicable, ensuring that the cases have been cleaned following the instructions above; where cases are not provided, package loosely in suitable pouch or cloth wrap. - Only legally marketed, and locally approved sterilization barriers (e.g., wraps, pouches or containers) should be used for packaging terminally sterilized devices, in compliance to the manufacturer's instructions. In the US, it is recommended that system trays be double wrapped using two FDA-cleared wraps prior to sterilisation. 	
Sterilisation	US	Outside US
	<ul style="list-style-type: none"> - Steam Sterilisation, Pre-vacuum Cycle - Temperature: 132°C (270°F) - Exposure time: 4 minutes - Dry-time: 40 minutes - Cool-time: 20 minutes 	<ul style="list-style-type: none"> - Steam Sterilisation, Pre-vacuum Cycle - Temperature: 134°C (273°F) - Exposure time: 3 minutes - Dry-time: 60 minutes - Cool-time: 30 minutes

Other sterilisation methods are possible but must be validated beforehand. Automatic cleaners and autoclaves must be validated by the hospital and regularly checked to guarantee the recommended sterilisation temperatures are reached for the entire exposure time.

If sterilisation containers with paper filters are used, it is advisable to use a new filter for each sterilisation. If after having followed this sterilisation method there is still water in the sterilisation containers or on/inside the device, the device must be dried and sterilisation repeated.

WARRANTY STATEMENT

Enztec Limited devices are manufactured for use only by qualified medical personnel who are trained in their use. All Enztec Limited devices are warranted to be free from defects in workmanship and materials for one (1) year from the date of sale. Any Enztec Limited device with a defect during the applicable warranty period will be repaired or replaced. Enztec Limited shall not be liable, expressly or impliedly, for:

- a. Any damages which arise or are caused, whether by the customer or by any of the users of the device or equipment, as a result of
 - (i) misuse, mishandling, and/or improper operation
 - (ii) repairs, modifications, or alterations performed by any person or entity other than Enztec Limited, or their authorised representatives
 - (iii) incorrect or incomplete inspection, cleaning and/or maintenance, or
 - (iv) use in combination with adaptors and/or equipment, or use in any manner or medical procedure, other than those for which it is designed; and
- b. Any special, indirect, and/or consequential damages of any kind and however caused arising from the sale or use of the device and equipment.

THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, AND/OR STATUTORY, INCLUDING, BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY, FITNESS AND/OR SUITABILITY FOR A PARTICULAR PURPOSE, AND ALL OTHER OBLIGATIONS OR LIABILITIES ON "ENZTEC LIMITED" BEHALF.

Enztec Limited neither assumes nor authorises any person to assume for it any other liabilities in connection with the sale of said devices and equipment. To ensure proper use, handling, and care of devices and equipment, consult the applicable catalogue, brochure, instruction manual, teaching film, and other literature which is included with the product and/or otherwise available from the company, upon request.

For further information relating to the use of this device or complaints please contact your Enztec Limited representative or distributor.

Rx: Federal law restricts this device to sale or use by or on the order of a physician.

SYMBOL TRANSLATION - see enztec.com/ifu	
QTY	Quantity
COUNTRY OF ORIGIN	Country of Origin

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Advena Ltd.