


# Classic Patella Saw Guide

An elegantly simple resection guide.

With a light-weight all-in-one design, the Enztec Classic Patella Saw Guide is designed for efficiency and flexibility.

The in-built stylus gives you the flexibility to choose from three resection heights. The efficient locking mechanism allows for rapid clamping while providing you with confident and consistent resections.

You can rely on the Enztec Classic Patella Saw Guide for easy and accurate resection.



## Features

- Guides a 1.27mm saw blade
- Multi-level stylus for resection height (8mm, 9mm, 10mm)
- Ratcheting locking mechanism
- Wide excursion, dual sided saw slot
- Easy to clean with no removable parts

## Specifications

- Materials: Stainless Steel
- Weight: 370g / 0.8lbs
- Length: 230mm / 9.0in
- Add your logo: available