



# Femoral Head Pusher

Supporting a direct superior approach.

With limited surgical access, the Enztec Femoral Head Pusher will perform the dual function of a head impactor and head pusher.

The longer, thinner impaction shaft allows better control and pressure to be maintained on the Femoral Head during movement and provides enhanced visibility.



## Features

- Compatible with 28, 32, 36 and 40mm femoral heads
- Replaceable head
- Thin design for increased visibility
- Easy to clean removable head that can be quickly reconnected

## Specifications

- Materials: Silicone, Stainless Steel, Propylux
- Weight: 200g / 0.4lbs
- Handle Length: 244mm/9.6in