




Femoral Impactor

Providing full control of femoral component placement and impaction.

With over 15 years of clinical history, the Enztec Femoral Impactor is easy to use and can be configured to work with any implant system. The robust design securely holds your entire size range of femoral components, while minimising interference with bone and soft tissue.

With over 18,000 units in the market, the Enztec Femoral Impactor is one of our most recognised and loved instruments. The Enztec Femoral Impactor is already compatible with a wide range of femoral components.



Compatible implants

- ZB Vanguard
- Corin Unity & Apex
- Microport Evolution
- Exactech Truliant
- + more

Features

- Intuitive to use
- Configurable to suit your femoral components
- Easy to clean
- Replaceable impaction head

Specifications

- Materials: Stainless Steel, Propylux
- Weight: 630g / 1.4lbs
- Length: 225mm / 8.8in
- Add your logo: available