




Classic Pin Puller

An elegantly simple solution for pin extraction.

The Enztec Classic Pin Puller provides simple and reliable Steinmann pin extraction. The innovative labyrinth design allows pins with minimal protrusion to be securely grasped, giving you confidence to extract pins in the most difficult circumstances.

With a novel, single-piece design, the Enztec Classic Puller is simple to use and easy to clean.

The Enztec Classic Pin Puller is a lightweight, reliable pin puller you can count on.



Features

- Compatible with 3.2mm / 0.125in Steinmann pins
- Grasps pins with only 3.2mm / 0.125in exposed
- Contoured handle provides comfort and grip
- Patented design (US D796,039)

Specifications

- Materials: Stainless Steel
- Weight: 120g / 0.3lbs
- Length: 145mm / 5.7in
- Add your logo: available