




Classic Pin Puller with Anvil

Providing additional support for pin extraction.

The Enztec Classic Pin Puller with Anvil provides simple and reliable Steinmann pin extraction. The integrated anvil allows pins to be impacted and provides enhanced grip when extracting stubborn pins.

The innovative labyrinth design allows pins with minimal protrusion to be securely grasped, giving you confidence to extract pins in the most difficult circumstances.

With a novel, single piece design, you can rely on the Enztec Classic Pin Puller with Anvil for simple pin impaction and extraction.



Features

- Compatible with 3.2mm / 0.125in Steinmann pins
- Grasps pins with only 3.2mm / 0.125in exposed
- Ergonomic handle design
- Integrated impaction anvil
- Patented design (US D796,039)

Specifications

- Materials: Stainless Steel
- Weight: 140g / 0.3lbs
- Length: 145mm / 5.7in
- Add your logo: available